GE Healthcare

Ettan IPGphor 3

Safety Handbook
Sicherheitshandbuch
Manual de seguridad
Manuel de sécurité
Manuale di sicurezza
Säkerhetshandbok

Ettan

Important user information

All users must read this entire handbook to fully understand the safe use of Ettan IPGphor 3.

WARNING!

The following **WARNING** signs highlight instructions that must be followed to **avoid personal injury**. It is important not to continue until all stated conditions are met and clearly understood.

WARNING! Lethal high voltage above 1 kV.

WARNING! Read the instructions to avoid hazardous conditions.

CAUTION!

The **Caution** sign highlights instructions that must be followed to **avoid damage to the product or other equipment**. It is important not to continue until all stated conditions are met and clearly understood.

Note

The **Note** sign is used to indicate information important for trouble-free and optimal use of the product.

CE certification

This product meets the requirements of applicable CE-directives. A copy of the corresponding Declaration of Conformity is available on request.

The **CE** symbol and corresponding declaration of conformity is valid for the instrument when it is:

- used as a stand-alone unit, or
- connected to other CE-marked GE Healthcare instruments, or
- connected to other products recommended or described in this handbook, and
- used in the same state as it was delivered from GE Healthcare except for alterations described in this handbook.

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1 Safety information

1.1 Introduction

IMPORTANT! All users must read this entire handbook to fully understand the safe use of Ettan™ IPGphor™ 3.

IMPORTANT! The Ettan IPGphor 3 system is intended for research use only and should not be used in any clinical or in Vitro procedures for diagnostic purposes.

IMPORTANT! To avoid any risk of injury, operation and user maintenance must be performed by appropriately trained personnel only, and in accordance with the instructions.

The purpose of this Safety Handbook is to describe:

- Safety precautions
- Emergency procedures
- Recycling of materials

1.1.1 Ettan IPGphor 3

The Ettan IPGphor 3 Isoelectric Focusing System performs first dimension isoelectric focusing (IEF) in 2-dimensional (2-D) electrophoresis.

The system comprises:

- The Ettan IPGphor 3 unit, combined with separate strip holders or a manifold
- Ettan IPGphor 3 Control Software
- Immobiline™ DryStrip gel strips, which contain an immobilized pH gradient (IPG)

Fig 1-1. The Ettan IPGphor 3 unit.

1.1.2 Associated documentation

- The *Ettan IPGphor 3 Instrument User Manual* describes the operation of Ettan IPGphor 3. It contains technical information, instructions for troubleshooting and user maintenance.
- The *Ettan IPGphor 3 Software User Manual* contains instructions for controlling Ettan IPGphor 3 Isoelectric Focusing Unit from a PC.
- The *Ettan IPGphor Cup Loading Manifold User Manual* contains instructions for first-dimension isoelectric focusing of proteins on IPG strips.
- Instructions attached to each package of Immobiline DryStrip.

1.1.3 Warning symbols

WARNING! The Warning symbol highlights instructions that must be strictly followed to avoid personal injury. Be sure not to proceed until the instructions are clearly understood and all stated conditions are met.

1.1.4 Caution notices

CAUTION! The Caution sign highlights instructions that must be followed to avoid damage to the product or other equipment. Be sure not to proceed until the instructions are clearly understood and all stated conditions are met.

1.1.5 Note signs

Note: *Indicates information important for trouble-free and optimal use of the product.*

1.2 *Safety precautions*

- Read the user manuals before using the Ettan IPGphor 3.
- The system is designed for indoor use only.
- Any equipment connected to the system should meet the requirements of the EN 61010-1 or other international safety standard. Within the EU, connected equipment must be CE-labeled.
- Select a place for Ettan IPGphor 3 where it is possible to access the mains power switch on the rear panel of the instrument.
- Do not operate the Ettan IPGphor 3 instrument in any way other than that described in the Ettan IPGphor 3 user manuals.
- For cleaning instructions, refer to *Ettan IPGphor 3 Instrument User Manual*.
- For installation of the PC, refer to the manufacturer's manual.

1.3 Warnings

1.3.1 Instrument operation

WARNING! Ettan IPGphor 3 is a high voltage instrument that can cause fatal electrical shock if the safety features are disabled. The safety lid must be securely closed before starting a protocol.

WARNING! The instrument covers must not be opened by the user. It contains electrical circuits which can give a lethal electric shock. Service and planned maintenance should be performed by personnel authorized by GE Healthcare.

WARNING! Plug the instrument into a properly grounded outlet.

WARNING! The vents at the front, back, and bottom of the instrument must not be obstructed. Air must be free to circulate for the cooling apparatus to function correctly. There should be a minimum of 10 cm unobstructed space behind the rear panel of the instrument.

WARNING! Always disconnect the power cord before servicing.

WARNING! This is a Class A product. In a domestic environment, it might cause radio interference, in which case the user might be required to take appropriate measures.

WARNING! Avoid spilling liquids on the body of the instrument. If large volumes of liquid have penetrated the casing of the instrument and come into contact with the electrical components, immediately switch off the instrument and contact an authorized service technician.

1.4 Cautions

1.4.1 Instrument operation

CAUTION! Exceeding the recommended current limit of 50 μ A per IPG strip can cause the strip to burn and may damage the instrument.

CAUTION! During isoelectric focusing, do not lean on the safety lid, do not apply excess pressure or uneven weight to the lid, and do not place any items on the lid. Such pressure could cause arcing between the strip holder electrodes and the electrode areas, damaging the instrument.

CAUTION! Operate Ettan IPGphor 3 on a level surface.

CAUTION! Avoid scratching the electrode areas.

CAUTION! Only use strip holders, manifolds, and spare parts supplied or specified by GE Healthcare.

1.4.2 IPG strips, strip holders, and reagents

CAUTION! The ceramic strip holders and manifold trays are very brittle and fragile. Handle carefully: Do not bump or drop.

CAUTION! Always wear gloves when handling IPG strips and the equipment that comes in contact with them. This will help minimize protein contamination which can result in artifactual spots in stained 2-D gel patterns.

CAUTION! Clean strip holders and manifold with the strip holder cleaning solution provided. Clean all other components that come in contact with the IPG strip or the sample with a detergent designed for glassware. Rinse well with distilled water.

CAUTION! Use the appropriate rehydration volume for the IPG strip length.


CAUTION! Solutions containing urea may be warmed briefly. Do not heat any solutions containing urea above 30–40 °C as isocyanate, a urea degradation product, will carbamylate the proteins in the sample, thus changing their isoelectric points.

CAUTION! All chemicals should be of the highest purity. Double distilled water should be used.

1.5 Labels and indicators

1.5.1 Rating label

The rating label (see figure below) is located on the instrument's rear panel.

1.5.2 Safety label

Safety labels (see figure below) are attached to both the light protective cover and the transparent lid of the safety lid. The safety label warns the user of risk for personal injury. Do not proceed until the instructions are clearly understood and all stated conditions are met.

1.5.3 Indicator and switch on the instrument unit

The Ettan IPGphor 3 instrument is equipped with the following indicator and switch:

Indicator/Switch	Color	Description
HV ON (indicator lamp)	Green	The HV ON lamp indicates that the high voltage power supply is on and isoelectric focusing is in progress. This lamp is not illuminated during the rehydration phase of a protocol.
Power (switch)	–	Switches on/off power to the instrument. Located on the rear panel, see Fig 1-2.

1.6 Emergency procedures

1.6.1 Emergency shutdown

In a situation where there is a risk of injury, turn off the mains power by switching the mains power switch (located on the rear panel) to position 0.


Fig 1-2. Mains power switch on the rear panel.

1.6.2 Power failure routine

In the event of a power failure, the run is interrupted in an undefined state.

1.6.3 Restart procedure

In the event of system shutdown due to power failure, emergency stop or process interruption, malfunctions must be rectified before the Ettan IPGphor 3 is restarted.

To restart Ettan IPGphor 3, follow the start-up instructions in *Ettan IPGphor 3 Instrument User Manual*.

1.6.4 Replacing fuses

Refer to instructions in *Ettan IPGphor 3 Instrument User Manual*.

1.7 Recycling

Ettan IPGphor 3 contains the following materials:

- Stainless steel
- Aluminum
- Platinum
- Titanium
- Ceramic
- Plastic:
 - Acrylic
 - Polycarbonate (PC)
 - Polyethylene (PE)
 - Polymethylmetakrylat (PMMA)
 - Polyurethane (PUR)
- Silicone
- Electronic components

When taking Ettan IPGphor 3 out of service these different materials must be separated and recycled according to local regulations, see also section 1.7.1 below.

1.7.1 Waste Electrical and Electronic Equipment (WEEE)

This symbol indicates that the waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Please contact an authorized representative of the manufacturer for information concerning the decommissioning of your equipment.

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www.gehealthcare.com

GE Healthcare Bio-Sciences AB
Björkgatan 30
751 84 Uppsala
Sweden

Ettan, Immobiline, and IPGphor are trademarks of GE Healthcare Ltd. GE tagline and GE monogram are trademarks of General Electric Company.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. GE Healthcare reserves the right, subject to any regulatory and contractual approval, if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your local GE Healthcare representative for the most current information.

© 2005 General Electric Company – All rights reserved.

GE Healthcare Bio-Sciences AB, a General Electric company.

GE Healthcare
Munzinger Strasse 9, D-79111 Freiburg, Germany

GE Healthcare
Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare
800 Centennial Avenue, P.O. Box 1327, Piscataway, NJ 08855-1327, USA

GE Healthcare
Sanken Bldg. 3-25-1, Hyakunincho, Shinjuku-ku, Tokyo 169-0073, Japan

Asia Paci.c Tel: +852 2811 8693 Fax: +852 2811 5251 • Australasia Tel: + 61 2 9899 0999 Fax: +61 2 9899 7511 • Austria Tel: 01/57606-1619 Fax: 01/57606-1627 • Belgium Tel: 0800 73 888
Fax: 03 272 1637 • Canada Tel: 800 463 5800 Fax: 800 567 1008 • Central, East, & South East Europe Tel: +43 1 982 3826 Fax: +43 1 985 8327 • Denmark Tel: 45 16 2400 Fax: 45 16 2424 • Finland &
Baltics Tel: +358 (0)9-512 39 40 Fax: +358 (0)9 512 39 439 • France Tel: 01 69 35 67 00 Fax: 01 69 41 96 77 • Germany Tel: 0761/4903-490 Fax: 0761/4903-405 • Italy Tel: 02 27322 1 Fax: 02 27302 212
• Japan Tel: +81 3 5331 9336 Fax: +81 3 5331 9370 • Latin America Tel: +55 11 3933 7300 Fax: +55 11 3933 7304 • Middle East & Africa Tel: +30 210 9600 687 Fax: +30 210 9600 693 • Netherlands
Tel: 0165 580 410 Fax: 0165 580 401 • Norway Tel: 815 65 555 Fax: 815 65 666 • Portugal Tel: 21 417 7035 Fax: 21 417 3184 • Russia & other C.I.S. & N.I.S Tel: +7 (095) 232 0250, 956 1137 Fax: +7 (095)
230 6377 • South East Asia Tel: 60 3 8024 2080 Fax: 60 3 8024 2090 • Spain Tel: 93 594 49 50 Fax: 93 594 49 55 • Sweden Tel: 018 612 1900 Fax: 018 612 1910 • Switzerland Tel: 0848 8028 12 Fax:
0848 8028 13 • UK Tel: 0800 616928 Fax: 0800 616927 • USA Tel: 800 526 3593 Fax: 877 295 8102

imagination at work