OneTouch Verio Pro+

Johnson & Johnson
MEDICAL

ONETOUCH[®]
VerioPro[™]+

Cordis

● DePuy

ETHICON

 ETHICON ENDO-SURGERY

LIFESCAN

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Performance of a New Glucose Meter Designed for Minimizing Infection Risk

L. Jones¹, G. Teft², J. J. Mahoney¹, J. M. Ellison¹. ¹LifeScan, Inc., Milpitas, CA, ²LifeScan Scotland, Inverness, United Kingdom

According to the Centers for Disease Control and Prevention (CDC), finger-stick devices, blood glucose meters, and the hands of health care workers are potential vehicles for the indirect transmission of blood borne pathogens in clinical settings. Advancements in technology for Assisted Monitoring of Blood Glucose (AMBG) might reduce the risk of bloodborne pathogen transmission among persons receiving diabetes care. The OneTouch Verio™Pro+ glucose meter incorporates an onboard test strip ejector that facilitates hands-free test strip disposal and minimizes health care provider (HCP) contact with blood-dosed test strips. This feature, along with the meter's compatibility with common disinfectants, may enhance infection control practices and may help reduce the risk of infection in clinical care settings.

In addition to infection control, glucose meters used for AMBG also require quality control (QC) tests to be performed regularly. The OneTouch Verio™Pro+ meter provides on-screen prompts to remind HCPs to perform QC tests and automatically recognizes when control solution is applied. After reading the Operator's Guide, 100% (38/38) of assessed HCPs reported that they knew how to correctly respond to a control solution prompt and 89.5% (34/38) stated that reminders to conduct a QC test will help them to be compliant with their institution's regulatory requirements.

Table 1 shows the results of a study conducted to ascertain meter performance at 3 clinics and 1 hospital with 189 capillary samples (hematocrit range: 22.9-52.1%), 177 venous samples (hematocrit range: 22.9-54.7%), and 200 arterial samples (hematocrit range: 23.0-59.8%) compared to YSI reference values. The OneTouch Verio™Pro+ glucose meter shows excellent analytical performance across the glucose range and incorporates many important clinical attributes required for AMBG and clinical settings.

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FAD-GDH and other chemical reagents
• Dried onto the Pd layer in the sensing chamber.

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Endogenous compounds

Interferent	High Normal Endogenous Concentration ² (mg/dL)	High Normal Endogenous Concentration ² (mmol/L)	Maximum Allowed Concentration ¹ (mg/dL)	Maximum Allowed Concentration ¹ (mmol/L)	Pass/Fail
Ascorbic Acid	2.0	0.11	6.0	0.34	Pass
Bilirubin	1.20	0.021	20.0	0.34	Pass
Glutathione	0.535	0.017	68.4	2.22	Pass
Uric Acid	8.0	0.48	9.97	0.59	Pass
Cholesterol	200	5.17	700	18.1	Pass
Creatinine	1.3	0.11	5.1	0.45	Pass
Haemoglobin	200	0.031	200	0.031	Pass
Lactic Acid	19.8	2.20	59.4	6.59	Pass
*Triglycerides	250	2.82	3000	33.9	Pass
Urea	85.8	14.3	257.4	42.9	Pass

*Endogenous substances are produced or originated from within the body.

‡High normal concentration is the level of endogenous substance observed in majority of patients.

†NO significant interference to the accuracy of the OneTouch® Verio™ Pro+ System was found at the concentration shown

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*Exogenous substances originate from outside the body, such as medications or vitamins.

‡High therapeutic concentration is the level at which that drug is expected to be effective without causing any serious problems to the patient.

† NO significant interference to the accuracy of the OneTouch® Verio™ Pro+ System was found at the concentration shown

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Exogenous compounds

Interferent	High Therapeutic Exogenous Concentration ² (mg/dL)	High Therapeutic Exogenous Concentration ² (mmol/L)	Maximum Allowed Concentration ¹ (mg/dL)	Maximum Allowed Concentration ¹ (mmol/L)	Pass/Fail
Mannitol	600	32.9	1800	98.8	Pass
Metformin	0.40	0.031	4.0	0.31	Pass
Methyl Dopa	0.75	0.036	1.5	0.071	Pass
Phenytoin	2.0	0.079	5.0	0.20	Pass
Pioglitazone	1.5	0.042	4.5	0.13	Pass
Procainamide	1.20	0.051	2.40	0.10	Pass
Quinidine	0.60	0.019	1.20	0.037	Pass
Salicylate	30	1.87	60	3.75	Pass
Simvastatin	2.7	0.065	8.1	0.19	Pass
Tetracycline	0.5	0.011	1.51	0.034	Pass
Tolazamide	3.31	0.11	9.93	0.32	Pass
Tolbutamide	10.8	0.40	64	2.37	Pass
Warfarin	0.31	0.010	1.0	0.032	Pass
Galactose	5.0	0.28	15	0.83	Pass
Lactose	0.5	0.014	4.1	0.11	Pass
Maltose	120	3.33	360	9.99	Pass
Sucrose	2670	78.0	5430	158.6	Pass

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+ THE NEW ONETOUCH[®] VERIO[™] PRO+

Because professionals
deserve a professional system

OneTouch® Verio™Pro+ System Quick Start Guide

Provides health care professionals in a multi-patient setting with safe, easy, accurate testing

Perform a blood glucose test*

1. Insert test strip to turn meter on

Insert two prongs of a test strip into test strip port with gold side and silver prongs facing you.

No coding required.

Apply Blood screen appears.

2. Get a blood sample*

Clean the puncture site with clean, soapy water, rinse, and dry.

Puncture finger with a single-use lancing device and gently squeeze the finger until a round drop of whole blood appears.

3. Apply to test strip (fingertip or syringe)

Apply blood sample from fingertip to either side of test strip.

Remove the metal needle from the syringe. Apply blood sample from a syringe to either side of test strip.

4. Read result

Meter displays result.

When **View Previous** appears, you can:

- > Press to see previous result
- > Then press to hide previous result

5. Remove used test strip

Point meter down into a biohazard container and press to release the test strip.

Disinfect meter after each patient

Step 1

Use a soft cloth dampened with the recommended disinfectant solution† or a presaturated wipe to completely pre-clean the meter surface.

Step 2

Use a new cloth, prepared as in Step 1. Squeeze out excess liquid before dampening outer meter surfaces. Ensure that the surface remains damp for the prescribed time as per the disinfectant manual.

*For venous or arterial samples, refer to the OneTouch® Verio™Pro+ Operator's Guide.

†Ten percent bleach solution or wipes containing a minimum 0.55% sodium hypochlorite. For more information, refer to OneTouch® Verio™Pro+ Operator's Guide.

ONETOUCH® Verio™Pro+

Review the OneTouch® Verio™Pro+ Operator's Guide before using the meter.

OneTouch® Verio™Pro+ System Quick Start Guide

Provides health care professionals in a multi-patient setting with safe, easy, accurate testing

Perform a control solution test

1. Identify level

Control solution level is printed on each label and indicated by the cap and label color.

Do not use after the expiration or discard date, whichever comes first.

2. Insert test strip and prepare the control solution

Insert two prongs of a test strip into test strip port with gold side and silver prongs facing you.

No coding required.

Apply Blood appears. Discard the first drop of the control solution.

3. Apply sample to strip

Squeeze the second drop of control solution into well on cap top. Touch test strip channel to solution until channel fills completely. The meter automatically detects the control solution.

4. Select control solution level

Meter counts down from 5 to 1 and displays three levels. Press **▲** or **▼** to select the level used and press **OK**.

If **Please Verify** screen appears:
> You selected the wrong level. Select **Change Level**.

OR

> The result does not match the control solution level you used. Select **Confirm** and test again. If the screen appears again, **do not** use the meter, test strips, or control solution.*

5. Read result

Once a control solution level is selected, the meter will display the control solution result along with the level of control solution selected.

Press **▲** to release the test strip.

Set control solution reminders

Step 1

Press and hold **OK** to turn the meter on, then press **OK** to go to main menu, then select **Settings**.

Highlight **Meter Care** and press **OK**.

Step 2

Select **Control Reminder** and press **OK**.

Step 3

Press **OK** with **Reminder** highlighted.

Step 4

Select **On** and press **OK**.

Step 5

Highlight **No. of Days** and press **OK**, then press **▲** or **▼** to select the desired number of days in between each control solution test and press **OK** again.

Step 6

Highlight **Time** and press **OK**, then press **▲** or **▼** to select the desired time the reminder will appear and press **OK** again.

Reminder

Within an 8-hour window after the desired reminder settings, a **Reminder** appears when the meter is turned on.

*If Please Verify screen continues to appear, contact Customer Service or a LifeScan sales representative.

ONETOUCH® VerioPro+

Thank you